

In-gel Digestion protocol for Mass Spec Analysis

Wash two sets of 500 µl snap cap micro centrifuge tubes prior to use with 2X 500 µl Methanol and 500 µl Milli-Q water rinse. Cut the protein band as small as possible using razor and take the band into the micro centrifuge tube.

A. RINSING

- 1) If gel slices are still wet with destain /etc. Rinse with 25 mM ammonium bicarbonate (ABC).

B. WASHING AND DEHYDRATION

- 1) Dehydrate gel with 50 µl of solution [A] ([A] = 2:1 mixture of Acetonitrile (ACN): 50 mM ABC) for 5 min. ABC = Ammonium Bicarbonate
- 2) Remove supernatant and add 50 µL of 25 mM ABC for 2 min.
- 3) Dehydrate gel slice with 50 µL of solution [A] for 5 min.
- 4) Remove supernatant and add 50 µl of 25 mM ABC for 2 min.
- 5) Dehydrate gel slice with 50 µL of solution [A] for 5min.

C. REDUCTION/ALKYLATION

- 1) After step B (Two or three cycles of dehydration/rehydration) rehydrate gel slice in 20-50 µL of 10-20 mM Dithiothreitol (DTT)
- 2) Place gel slice at 60 °C for 1 hour
- 3) Remove DTT, and rinse gel slice with 50 µl 25 mM ABC
- 4) Add 20-50 µL of 50-100 mM Iodoacetamide (IAA).
- 5) Incubate gel slice at RT in the dark for 20 min.
- 6) Remove IAA, and rinse gel slice with 50 µL ABC.

D. DEHYDRATION (Continued)

- 1) Dehydrate gel slice with 50 µL solution A for 5 min.
- 2) Remove supernatant and add 50 µl of 25 mM ABC for 2 min.
- 3) Optional (repeat dehydration and rehydration 1 more time)
- 4) Remove final supernatant and dry spots in a speed vac (ca. 5-15min).

E. DIGESTION

- 1) Prepare trypsin stock by adding 1 mL of ice-cold 1 mM HCl to a standard 20 µg trypsin vial (Sigma/Promega MS grade trypsin). Store on ice. Make aliquots of 5 µL. For each gel band, add 95 µL of 25 mM ABC to the stock (concentration is 100 ng)
- 2) Rehydrate the gel slice with trypsin, (100 µL) at 0 °C on ice.
- 3) Incubate gel slice on ice 20-30 min until trypsin is absorbed.
- 4) After gel slice is completely rehydrated, add just enough 25 mM ABC to cover gel slice in tube
- 5) Incubate at 37 °C for at least 16 hours, typically overnight. (24 hrs recommended)
- 6) Periodically check gel slice to make sure it is still just covered with 25 mM ABC, add more if needed.

F. EXTRACTION

- 1) Remove supernatant, and vortex gel slice with 2x extractions of 100 μ l extraction solution (50:50 ACN:H₂O, 0.1% Trifluoroacetic Acid), add these extractions to supernatant. Collect supernatant of step 1, 2 and 3 and reduce the volume in speed vac. (Total volume approximately 500 μ l)
- 2) Speed vac the mixture to complete dryness and resuspend in 10 μ L of 0.1% TFA in water. This solution is used for mass spec after ZIP Tip.

REAGENTS PREPARATION:

- 1) Preparation of 25 mM ABC
Add 0.1977g NH₄HCO₃ to 100 ml Milli-Q water and mix until dissolved.
- 2) Preparation of 50 mM ABC
Add 0.395g NH₄HCO₃ to 100 ml Milli-Q water and mix until dissolved.
- 3) Preparation of 10 mM DTT
Add 1.5 mg DTT to a 1.5 ml Eppendorf tube.
Add 1 ml of 50 mM ammonium bicarbonate and mix until dissolved.
- 4) Preparation of 50 mM Iodoacetamide
Add 10 mg IAA to a 1.5 ml Eppendorf tube.
Add 1 ml of 50 mM ammonium bicarbonate.
- 5) Preparation of 0.1% TFA
Add 100 μ l of TFA to 99.9 ml Milli-Q water. Mix it.

CHEMICALS AND CONSUMABLES:

- 1) Ammonium Bicarbonate: Sigma Cat # A6141
- 2) Trifluoroacetic acid: Sigma Cat # 299537-100G
- 3) Acetonitrile: Sigma Cat # 34967
- 4) Trypsin: Sigma Cat # T6567
- 5) Zip tips: Sigma Cat # Z720070-96EA
- 6) IAA: Sigma Cat # I6125-25G
- 7) DTT: Sigma Cat # D9779-5G
- 8) Promega Trypsin: Promega Cat # V5111

Zip Tip C18 Protocol for MS Analysis:

I. Materials

- **Wetting solution** (100% HPLC-grade Acetonitrile)
- **Washing solution** (0.1% TFA in Milli Q water)
- **Sample preparation:** Adjust sample to 0.1% TFA
- **Elution solution** (50% Acetonitrile in 0.1% TFA):
300 μ L HPLC-grade Acetonitrile, 240 μ L Milli Q water, 60 μ L 1% TFA in Milli Q water

II. Procedure

Note: Resin bed provides back pressure, so set pipette to 10 μ L, depress plunger to dead stop and slowly release or dispense plunger throughout operation.

1. Equilibrate:

Aspirate 10 μ L wetting solution into tip and dispense to waste. Repeat. Aspirate washing solution into tip and dispense to waste. Repeat.

2. Bind & wash:

Bind peptides to Zip Tip pipette tip by aspirating and dispensing 8-10 cycles (simple mixtures), up to 15 cycles (complex). Aspirate **washing solution** and dispense to waste. Repeat wash 5 times.

Note: A 5% methanol in 0.1% TFA/water wash can improve desalting efficiency.

3. Elute:

Dispense 10 μ L of **elution solution** into clean 0.5 mL Eppendorf micro centrifuge tube using a standard pipette tip. Aspirate and dispense eluent through Zip Tip at least 8-10 times without introducing air.

Reference:

Shevchenko, A.; Tomas, H.; Havlis, J.; Oslen, J. V.; Mann, M. In-gel digestion for mass spectrometric characterization of proteins and proteomes. *Nature Protocols* (2006) Vol-6, 2856-2860.