

## Digestion protocol for single protein

This protocol for digestion of single protein the procedure follows.

- 1) Take 100  $\mu\text{g}$  of the individual protein.
- 2) Add 100  $\mu\text{l}$  of (100mM) Ammonium bicarbonate (ABC).
- 3) Add 20% of ACN (approx. 20  $\mu\text{l}$  of ACN).
- 4) Incubate for 45 minutes at 60°C.
- 5) **Reduction:** Add 5  $\mu\text{l}$  of (50mM) Dithiothreitol (DTT) and incubate for 45 minutes at 60°C (If protein do not have disulphide bonds, skip this step)
- 6) **Alkylation:** Add 5  $\mu\text{l}$  of (50mM) Iodoacetamide (IAA) and incubate in dark at 45 minutes at room temperature (If protein do not have disulphide bonds, skip this step)
- 7) Add Trypsin in 1:50 ratio (Enzyme:Substrate) and incubate for 24 hours at 37°C in shaking water bath.
- 8) Add 2  $\mu\text{l}$  Trifluoroacetic acid (TFA) or Formic acid to stop the Trypsin activity.
- 9) Then subject the sample to speed vacuum to complete dryness.
- 10) Add 10  $\mu\text{l}$  of 0.1 % TFA in water solution to reconstitute.
- 11) Then perform Zip-Tip to reduce the salt concentration from the sample.